PHOENIX MOLDED SHOES INSERTS & BRACES



PHOENIX

A New Erals Born

e at Phoenix Molded Shoes, Inserts, and Braces have built a company based on three principles. These principles are our foundation and are indispensable.

- 1) Great Customer Service (Our job is not complete until our product is dispensed and everyone involved is satisfied)
- 2) Quality Product
- 3) An Expeditious and Reliable Turn Around time

We work one on one with you, the practitioner, to provide you with the best fitting, most functional custom footwear and devices available today. We know that communication is the key, you know what you want for your patient and you can't always express that on an order form. That is why we are here to answer your calls or call you regarding orders or anything you would like to discuss.

Having owned and operated our own Orthotic, Prosthetic, and Pedorthic facilities we have been making and dispensing custom molded shoes, orthotics, and braces for over 25 years. This is why Padded Collars and Mild Rocker Soles always come standard on our shoes with the exception of some of our Premier Dress Styles; there is never an additional charge for these modifications because we know these are needed to help insure a proper fit.

We would like to thank you for putting your trust in our company this is not an issue we take lightly, we understand that there are a lot of choices out there for your custom needs and we will do everything in our power to provide you with the best quality product in the most timely manner possible. Once again thank you and we truly value your business.

Keith B. Brewer C.Ped., L.Ped.

Vertl & Breity

CEO





CUSTOM MODIFICATIONS



CLOSURES

Lace

Speed Lace

Hook Lace

Velcro Closures

- With or without D-rings
- Normal or Reverse

Zippers

Buckles

OPENINGS

Regular

Semi-Surgical (1" longer than regular)

Surgical (to toes)

LININGS

B-Foam

Leather

Plastazote®

SOLING

12 Iron Crepe

18 Iron Crepe

24 Iron Crepe

Heavy Duty Ribbed Sole

Wedge Lug

Lug w/Distinct Heel

Leather Sole with Distinct Heel

Leather Toe Tips

Full Leather Slip Sole (welted construction)

OUTSOLE MODIFICATIONS

Heel and Sole Flares

Heel and Sole Wedges

Wide Bases

Steel Shank (1/2 or full)

Caliper Installation (yours or ours)

Stirrup Installation (yours or ours)

Heel or Full Build-Ups (crepe or cork

with leather cover)

MISCELLANEOUS

Duplicate Cast (we hold cast for 3 months)

T-Strap

Reinforce Counters



INSOLES

1/4" Pink Plastazote® and 1/4" Silver Plastazote®

1/4" Pink Plastazote® and 1/8" Poron® and

1/4" Silver Plastazote®

1/2" Silver Plastazote®

1/2"Pink Plastazote®

1/2" Microcell Puff®

Spenco® Cover

Poron® Cover

Leather Cover

Toe or Forefoot Filler (steel shank required for trans-met to choparts)

Metatarsal Pads

Pockets on Insoles

Base Depressions

Reinforce for Brace



NEO-FLEX

Custom Molded Shoes

All of the qualities of a custom molded shoe with the advantages of a stretchable upper to aid in comfortably fitting feet with:

Hammer Toes
 Bunions
 Severe Edema
 HyperSensitivity

Neo-Flex can be substituted for leather in most of our styles.

The Neo-Flex option can be selected simply by placing an "N" after the style #.



The only company to offer stretchable uppers in 9 colors.



















Black

White

Navy

Brown

Cord

Luggage Tan

Taupe

Bone

COLORS & LEATHERS











White



Navy



Brown



Dk Brown



Cordovan



Luggage Tan



Taupe



Bone





Black



White



Taupe



Bone



Brown

Nивиск



Black



Brown



Taupe



Wheat

HEAVY DUTY LEATHER



Black



Brown



Cedar





CASUAL





Low Top Blucher Style Shoe with plain toe, single Velcro closure, and pillow back collar.



PH-1 Low Top Blucher Shoe with plain toe and lace closure



Low Top Blucher Shoe with plain toe and velcro closure





Low Top Blucher Style Shoe with decorative seamed toe and lace closure





CASUAL





CASUAL









High Top 5 1/2" Boot with plain toe and velcro closure (Also available in 4 1/2" Chukka PH-2CA)


with decorative seamed toe and lace closure (Also available in velcro closure PH-3CA)

CASUAL









PH-6

Low Top Blucher Style Shoe

with toe cap styling and lace closure (Also available in velcro closure PH-6A)

PH-16
Low Top Balmoral Style Shoe
with athletic styling and lace closure
(Also available in velcro closure PH-16A)



SPORT





Low Top Blucher Style Two Tone Sport Shoe with lace closure (Available in any combination of colors)



PH-25

Low Top Blucher Style Shoe with extended opening. Available in any combination of colors and leathers.



Low Top Blucher Style Shoe with extended opening and extra eyelets for added adjustability and the larger foot.



Low Top Blucher Style Golf Shoe with nipple soling and lace closure (Also available in velcro closure PH-18A)





Low Top Blucher Style Two Tone Sport Shoe with velcro closure (Available in any combination of colors)



Low Top Blucher Style Shoe with athletic styling and lace closure



High Top 5" Blucher Shoe with athletic styling and lace closure (Also available in velcro closure PH-4HA)



Low Top Blucher Style Shoe with athletic styling and velcro closure



PH-25A

Low Top Blucher Style Shoe with three Velcro straps and an extended opening.





WORK



Rugged 6" Workboot with double padded collar (Also available in velcro closure PH-7/6A)



Rugged 4 1/2" Chukka Height Workboot with padded collar (Also available in velcro closure PH-7CA)



Rugged 8" Workboot with double padded collar (Also available in velcro closure PH-7/8A)



Low Top Rugged Blucher Style
Work Shoe with plain toe and lace closure
(Also available in velcro closure PH-7LA)













Note: PH-11, PH-12 and PH-13 are Dress Styles and not recommended for neuropathic patients, patients prone to foot ulcerations, or patients who require a lot of room in the forefoot.



WOMEN'S DRESS



Mary Jane Style Women's Shoe with velcro closure



Women's Blucher Style Dress Shoe with patent leather accents (Available in black only)



Note: PH-9, PH-9B and PH-14 are Dress Styles and not recommended for neuropathic patients, patients prone to foot ulcerations, or patients who require a lot of room in the forefoot.

SANDALS












DIABETIC INSERT PACKAGES

 Pink Plastazote and Eva also available in tri-laminate with 1/16 poron between layers.









Trans-Met Filler





PPT Layer



Depression Modification











A Shoe is Born of



1. Receive cast.



3. Casts cure before modifying.



5. Vacuum form inserts.



7. Prepare cast for leather upper pattern.



2. Pour negative cast with plaster.



4. Shape and modify positive cast.



6. Grind and shape inserts.



8. Cut upper pattern.



Tradition and Craftsmanship





9. Stitch uppers.



11. Prepare for soling.



13. Finished shoes await inspection.



15. Polished.



10. Fit upper to cast.



12. Attach and shape soles.



14. Finished shoe checked for accuracy.



16. Shipped to you.



PHOENIX MOLDED SHOES INSERTS & BRACES

